

Q-Optics Through-the-Lens (TTL) Loupes

**NEW
ECLIPSE-TI FRAME**



Custom, uniquely designed and built to offer both comfort and precision performance.

Built with ultra-strong, lightweight materials, the Q-Optics TTL loupes are available in two styles: High Resolution and Micro Weight and are available in a variety of frames and colors.



HIGH RESOLUTION THROUGH-THE-LENS (TTL) LOUPES

Q-Optics Custom High Resolution TTL style Loupes feature an extra dense optical system. Multiple lenses and glass types combine to create up to a 25% increase in image clarity. Weight was not ignored however as the loupes still weigh on average 25 grams less than standard optical systems. TTL systems are custom made to align to the users facial features. High Resolution Through the Lens style loupes come in black with your choice of frame style and frame color.



MICRO WEIGHT THROUGH-THE-LENS (TTL) LOUPES

Q-Optics Micro Weight Through-the-Lens style Loupes feature a special lens spacing technique that uses both glass types and distance to help eliminate optical aberrations and bring an extreme reduction in weight. All Micro Weight TTL systems are custom aligned to the user and are available in a variety of frame styles and frame colors.

Custom, uniquely designed and built for comfort and precision performance.



CUSTOM TTL LOUPES

- Super light weight
- Four point ergonomic customization—no standard angles, widths, or loupe heights for a true custom fit
- Break-proof polycarbonate lenses
- Laser aligned convergent point
- True Rx—no prescription manipulation
- Titanium frame, safety frame, or sport frame
- Waterproof
- Military-grade anti-scratch coating
- Lifetime Guarantee on Titanium frames, optical seal, and bond of telescope to frame lens and cracking of frame lens
- 10" - 26" customized working distance

Available Frame Styles



Frame colors available:
 Pewter
 White
 Black
 Silver
 Cotton Candy
 Pearl Blue



Eclipse-ti



Frame colors available:
 Dark Pewter
 Pearl Pink
 Gunmetal Grey
 Copper



Safety Q3200



Frame colors available:
 Blue, Earth, Black



Standard Q867

Sizes: 54, 56

Frame colors available:
 Gun metal (Shown)

Fashion Q869

Sizes: 48, 50

Frame colors available:
 Rose (Shown)
 Black



Specifications

POWER	MODEL NUMBER	TYPE	WORKING DISTANCE	FIELD OF VIEW	COLOR	WEIGHT
Micro Weight						
2.5X	M-2.5 TTL	Galilean	10"-26" customized	Entire open mouth	Black	38.0g (1.3oz)
3.0X	M-3.0 TTL	Galilean	10"-26" customized	80-90% open mouth	Black	36.0g (1.3oz)
3.5X	M-3.5 TTL	Galilean	10"-26" customized	50-60% open mouth	Black	36.0g (1.3oz)
High Resolution						
2.5X	F-2.5 TTL	Galilean	10"-26" customized	Entire open mouth	Black	43.0g (1.5oz)

